# Operating manual


## VWR® IDP 761-C - Conductivity Sensor

EU cat. no

NA cat. no

663-0356

76470-834

CE

UK  
CA

**Legal Address of Manufacturer:**

Europe

VWR International BV  
Researchpark Haasrode 2020  
Geldenaaksebaan 464  
B-3001 Leuven  
+32 16 385011  
<http://be.vwr.com>

United States

VWR International LLC  
100 Matsonford Rd  
Radnor, PA 19087  
+1 800-932-5000  
[www.vwr.com](http://www.vwr.com)

Importer to UK:

VWR International Ltd  
Hunter Boulevard, Magna Park  
Lutterworth, Leicestershire, LE17 4XN  
<http://uk.vwr.com>


Country of origin: Germany.

EN Operating manual . . . . .	4
DE Bedienungsanleitung . . . . .	9
FR Mode d'emploi . . . . .	14
IT Manuale operativo . . . . .	19
PT Instruções de utilização . . . . .	24
ES Español . . . . .	29

## Overview

### Structure and function

#### Structure

#### Automatic sensor recognition

The sensor electronics with the stored sensor data is in the connecting head. The data include, among other things, the sensor type and series number. With each calibration, the calibration data is written in the sensor and the calibration history is recorded. The data is recalled by the meter when the sensor is connected and is used for measurement and for measured value documentation. Storing the calibration data in the sensor ensures that the correct cell constant is automatically used if the sensor is operated with several meters.

The digital transmission technique guarantees the failure-free communication with the meter even with long connection cables. If the sensor firmware is enhanced by VWR, it can be updated with the meter.

#### Recommended fields of application

- On-site measurements in rivers, lakes and waste water
- Fish farming
- Groundwater measurements
- Applications in the water lab

#### **NOTE**

**Concentrated or highly oxidizing acids and organic solvents can damage or destroy the sensor.**

## Measurement / Operation

### Commissioning

**Scope of delivery**

- Conductivity measuring cell IDP 761-C
- Operating manual

**Commissioning**

Connect the sensor to the meter.

### Cleaning

#### **NOTE**

**To clean the sensor, disconnect it from the instrument.**

**Outside cleaning**

We recommend to clean the sensor thoroughly, especially before measuring low conductivity values.

<b>Contamination</b>	<b>Cleaning procedure</b>
Lime sediments	Immerse in acetic acid for 5 minutes (volume share = 10 %)
Fat/oil	Clean with warm water that contains washing-up liquid

After cleaning, thoroughly rinse with deionized water and recalibrate if necessary.

**Aging of the conductivity measuring cell**

Normally, the conductivity measuring cell does not age. Special measuring media (e.g. strong acids and bases, organic solvents) or temperatures that are too high may considerably reduce its lifetime or lead to damage. The warranty does not cover failure caused by measuring conditions and mechanical damage.

**Disposal**

We recommend to dispose of the measuring cell as electronic waste.

## What to do if ...

Error symptom	Cause	Remedy
No temperature or conductivity display	<ul style="list-style-type: none"> <li>No connection between meter and conductivity measuring cell</li> <li>Cable defective</li> </ul>	<ul style="list-style-type: none"> <li>Establish connection between meter and conductivity measuring cell</li> </ul>
Measurement delivers implausible conductivity values	<ul style="list-style-type: none"> <li>Measuring range exceeded</li> <li>Contamination in the area of the electrodes</li> <li>Electrodes damaged</li> </ul>	<ul style="list-style-type: none"> <li>Make sure the correct sensor is being used for the application</li> <li>Clean the conductivity measuring cell (see section ).</li> <li>Return the sensor</li> </ul>
Incorrect temperature display	<ul style="list-style-type: none"> <li>The temperature sensor is not immersed deep enough in the measuring solution</li> <li>Temperature sensor defective</li> </ul>	<ul style="list-style-type: none"> <li>Observe the minimum immersion depth</li> <li>Return the conductivity measuring cell</li> </ul>


## Technical data

### General data

#### General features

Measuring principle	Four-electrodes measurement
Cell constant	$0.475 \text{ cm}^{-1} \pm 1.5 \%$
Temperature sensor	Integrated NTC 30 (30 k $\Omega$ at 25 °C / 77 °F)

#### Dimensions (in mm)

#### Weights

60 g (without cable)

#### Materials

Shaft	Epoxy
Connection head	POM
Conductivity electrodes	Graphite
Thermistor enclosure	Graphite

<b>Connection cable</b>	Length	1.5 m
	Diameter	4.3 mm
	Smallest allowed bend radius	Fixed installation: 20 mm Flexible use: 60 mm
	Plug type	Socket, 4 pins
<b>Pressure resistance</b>	Sensor with connection cable	IP 68 (2.5 x 10 <sup>5</sup> Pa or 2.5 bar)
<p>The IDP 761-C meets the requirements according to article 3(3) of the directive, 97/23/EC ("pressure equipment directive").</p>		
<b>Measurement conditions</b>	Conductivity measuring range	1 µS/cm ... 500 mS/cm
	Temperature range	-5 ... 70 °C (100 °C) 23 ... 158 °F (212 °F)
	Max. admissible overpressure	IDP 761-C: 2.5 x 10 <sup>5</sup> Pa (2.5 bar)
	Minimum depth of immersion	36 mm
	Maximum depth of immersion (at temperature)	Whole sensor + cable up to 70 °C (158 °F) Sensor shaft only (=120 mm) up to 100 °C (212 °F)
	Operating position	Any
<b>Storing conditions</b>	Recommended storing method	In air
	Storage temperature	0 ... 50 °C (32 ... 122 °F)
<b>Characteristics when delivered</b>	Temperature responding behavior	t <sub>99</sub> (99 % of the final value display after) < 20 s
	Accuracy of the temperature sensor	± 0.2 K

### Measuring ranges, resolution, accuracy

Measuring ranges, resolution	Measured parameter	Measuring range	Resolution
	æ [µS/cm]	0.0 ... 199.9	0.1
		200 ... 1999	1
	æ [mS/cm]	2.00 ... 19.99	0.01
		20.0 ... 199.9	0.1
		200 ... 500	1

Measured parameter	Measuring range	Resolution
$\rho$ (resistivity) [Ohm*cm]	2.00 ... 19.99	0.01
	20.0 ... 199.9	0.1
	200 ... 1999	1
$\rho$ (resistivity) [kOhm*cm]	2.00 ... 19.99	0.01
	20.0 ... 199.9	0.1
	200 ... 1999	1
$\rho$ (resistivity) [MOhm*cm]	2.00 ... 19.99	0.01
SAL	0.0 ... 70.0 according to the IOT table	0.1
TDS	0 ... 1999 mg/l	1
	2.00 ... 19.99 g/l	0.01
	20.0 ... 199.9 g/l	0.1
T [°C]	- 5,0 ... + 100,0	0.1

### Accuracy of the IDP measuring technique


Measured parameter	Accuracy ( $\pm 1$ digit)
$\rho$ , SAL, TDS	$\pm 0.5$ % of measured value
T [°C]	$\pm 0.1$

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VWR International AB
VWR International GmbH
VWR International, LLC