Operating manual


VWR® OPOX 11-3 - D.O. Sensor

EU cat. no

NA cat. no

664-0182

76470-824

CE

**UK
CA**

Legal Address of Manufacturer:

Europe

VWR International BV
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+32 16 385011
<http://be.vwr.com>

United States

VWR International LLC
100 Matsonford Rd
Radnor, PA 19087
+1 800-932-5000
www.vwr.com

Importer to UK:

VWR International Ltd
Hunter Boulevard, Magna Park
Lutterworth, Leicestershire, LE17 4XN
<http://uk.vwr.com>

Country of origin: Germany.

Contents


EN English	5
1 Overview	5
2 Measurement / Operation	6
3 Maintenance, cleaning, replacement	8
4 What to do if	11
5 Technical data	12
DE Deutsch	14
1 Überblick	14
2 Messen / Betrieb	15
3 Wartung, Reinigung, Ersatzbedarf	17
4 Was tun, wenn	21
5 Technische Daten	22
ES Español	24
1 Sumario	24
2 Medición / funcionamiento	25
3 Mantenimiento, limpieza, repuestos necesarios	27
4 Diagnóstico y corrección de fallas	31
5 Especificaciones técnicas	32
FR Français	34
1 Vue d'ensemble	34
2 Mesure / fonctionnement	35
3 Maintenance, nettoyage, pièces et fournitures de rechange	37
4 Que faire, si... ..	41
5 Caractéristiques techniques	42

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Structure

Automatic sensor recognition

The data of the sensor and sensor cap are recalled by the meter when the sensor is connected and is used for measurement and for measured value documentation. The calibration data are stored in the sensor. Thus, the calibration is automatically retained if the sensor is connected to another meter.

The digital transmission technique guarantees the failure-free communication with the meter even with long connection cables.

1.2 Recommended fields of application

Recommended fields of application

- On site measurements in rivers, lakes and wastewater
- Applications in water laboratories
- BOD measurements

2 Measurement / Operation

2.1 General information on the handling of the sensor cap

- Please do not touch the sensor membrane with your fingers unnecessarily
- Avoid any great mechanical stress of the sensor membrane (pressure, scratches).

2.2 Commissioning

Scope of delivery

- D.O. sensor OPOX 11-3 with sensor cap
- Check and storage beaker
- Operating manual

Getting the sensor ready for measuring

1.	Unscrew the protective hood from the sensor by hand.
2.	Pull off the protective cap (red).
3.	Put the protective hood on the sensor head and screw it tight by hand <u>as far as it will go</u> . The sensor is immediately ready to measure.

2.3 Measuring

Minimum immersion depth

Observe the required minimum immersion depth (see chapter 5 TECHNICAL DATA).

Approach flow

The OPOX 11-3 D.O. sensor enables precise measurements without any incident flow.

However, an incident flow of the sensor membrane improves the responding behavior of the sensor (see chapter 5 TECHNICAL DATA). The incident flow can be provided in different ways, e. g.:

- The flow of the water to be measured is sufficient (aeration tank, water pipe, stream)
- Slowly pull the sensor through the water by hand (lake, container), or
- Use a flow aid such as a magnetic stirrer

2.4 Function check and user calibration

Factory calibration

The OPOX 11-3 is factory calibrated. In the recommended application (see page 5), the measuring characteristics of the sensor cap remain stable for the specified service life. Thus, a user calibration is not usually required.


When does a function check or user calibration make sense?

A function check or user calibration can be useful in the following special cases:

- If the measured values appear to be implausible and it is assumed that the service life of the sensor cap is over
- Routinely within the framework of the company quality assurance

Check- or calibration medium

The user calibration takes place in water vapor-saturated air. The suitable conditions are easily provided with the aid of the calibration and storing beaker. For this purpose, moisten the sponge inside the calibration and storing beaker. Then insert the sensor in the calibration and storing beaker. The sensor membrane must be clean and dry for this.

Moisten the sponge in the calibration and storing beaker.

After inserting the sensor, wait for the temperatures of the sensor and storing cap to adjust.

Note

The steps of the user calibration are described in detail in the operating manual of the meter.

2.5 Storage

Always store the sensor with the calibration and storage vessel put on at a temperature in the range 0 ... + 50 °C (32 ... 122 °F).

3 Maintenance, cleaning, replacement

3.1 General maintenance instructions

Handling of the sensor cap

- Dirt and moisture under the sensor cap can affect the functioning and shorten the service life of the sensor cap. Therefore, make sure the working environment is clean and dry prior to removing the sensor cap.
- Please do not touch the sensor membrane unnecessarily with your fingers. Touch the sensor cap at the sides only.
- Avoid any great mechanical stress of the sensor membrane (pressure, scratches).
- Exposure to light, particularly daylight, of the interior of the sensor cap will, by-and-by, affect the measurement characteristics and shorten the service life of the sensor cap. Therefore, the interior of the sensor cap should not be exposed to direct sunlight. Avoid any exposure to light that exceeds the extent required for necessary maintenance and cleaning activities. Store dismantled sensor caps in a light-protected environment only.

3.2 Exchanging the sensor cap


Note

Before removing the sensor cap, observe the general maintenance instructions in section 3.1.


To exchange the sensor cap, proceed as follows:

Removing the sensor cap

1.	Pull the sensor out of the sample.
2.	Clean the outside of the sensor (see section 3.3.1).
3.	Unscrew the protective hood from the sensor by hand.
4.	Thoroughly clean and dry the sensor head once again.
5.	Grasp the sensor cap on the sides (arrows in figure on the left) and unscrew it from the sensor anti-clockwise.
6.	Remove the sealing ring. Caution: Do not push any tools or other sharp objects between the sealing surfaces. This would damage the sealing surfaces and sealing ring.
7.	Remove any remaining sealing grease from the seat of the sealing ring using a dry lens-cleaning towel. Wipe off any other deposits using a moist lens-cleaning towel. Wipe dry using a dry lens-cleaning towel.

Mounting the sensor cap

8.	Insert the sealing ring provided with the new sensor cap. Caution: Do not push any tools or other sharp objects between the sealing surfaces. This would damage the sealing surfaces and sealing ring. Do not touch the measuring window. Wipe any fingerprints off the measuring window as necessary.
9.	If necessary, clean the sealing ring once again.
10.	Cover the new sealing ring with a thin film of sealing grease. Remove any excess sealing grease using a lens-cleaning towel.
11.	Clean the surface of the measuring window using a lens-cleaning towel.
12.	Take the new sensor cap out of the container. Dry the hollow of the sensor cap completely using a soft cloth.
13.	Screw on the new sensor cap clockwise and by hand. The sealing ring is between the sensor cap and sensor. Caution: Do not use any tools. Do not screw on the sensor cap too tight.
14.	Put the protective hood on the sensor head and screw it tight by hand <u>as far as it will go</u> .
15.	Connect the sensor to the meter.
16.	Enter the cap coefficients (K1 to K5 and KC) on the meter. The cap coefficients are provided with the sensor cap (see operating manual of the meter).

3.3 Clean the sensor

3.3.1 Exterior cleaning

Extreme dirt on the sensor can affect the measuring characteristics. Biological deposits for example, consume oxygen and can, when occurring on the sensor membrane, impair the responding behavior and cause values that are too low. Therefore, we recommend regular visual inspections and exterior cleaning as necessary.

Pay attention to the following points for cleaning:

- First, thoroughly rinse the sensor with tapwater to remove loosely adhering dirt.
- Rough dirt on the sensor shaft can be brushed off with a soft brush. **Attention:** Do not use the brush in the area of the sensor membrane. Risk of damage!
- The sensor cap including the sensor membrane should be wiped with a soft and moist microfiber cloth.

- In the case of persisting dirt you can add some household washing-up liquid to the tapwater. Attention: Never use any alcohol for cleaning!

3.3.2 Interior cleaning of sensor cap and sensor head

If moisture or dirt have penetrated under the sensor cap, e.g. because the sensor cap is damaged, you can make the sensor ready for operation again as follows:

CAUTION

Only use **nonabrasive, alcohol-free** detergents, as otherwise the optical surfaces could be damaged.

1.	Remove the sensor cap (see section 3.2).
2.	Clean the sensor head and sensor cap: <ul style="list-style-type: none"> – Rinse all inner surfaces with tapwater – Remove contamination containing fat and oil with warm water and household washing-up liquid – Then rinse all inner surfaces with deionized water
3.	Pat dry all surfaces with a clean, lint free cloth.
4.	Allow the sensor and sensor cap to dry completely at a dry location so moisture can evaporate even from corners difficult to access. When doing, so protect the inside of the sensor cap from light.
5.	Put the sensor cap on (see section 3.2).

Note

If the sensor cap is visibly damaged it has to be replaced.

3.4 Checking the zero point of the sensor

There are two methods to check the zero point of the sensor:

- Measurement in a nitrogen atmosphere (recommended method)
- Measurement in a sodium sulfite solution according to DIN EN 25814/ISO 5814.

Test criterion

The sensor is OK if the measuring instrument displays < 0.5 % D.O. saturation after 15 minutes.

4 What to do if ...


Error symptom	Cause	Remedy
No temperature value or D.O. value	– No connection between meter and D.O. sensor	– Establish connection between meter and D.O. sensor
	– Cable defective	– Return the D.O. sensor
Measured value too high or too low or error message, <i>Error</i>	– Coating on sensor cap	– Clean the outside of the sensor (see section 3.3.1)
	– Membrane damaged	– Replace the sensor cap
	– Service life of the sensor cap over	– Check the sensor – Replace the sensor cap if necessary (see section 3.2)
	– Dirt inside the sensor cap and in the sensor head. – Protective hood not screwed on tight enough – Sensor cap untight or defective	– Dismantle the sensor cap – Clean the insides of the sensor cap and sensor head (see section 3.3.2) – Mount the sensor cap correctly and tighten the protective hood as far as it will go (see section 3.2) – Replace a defective sensor cap if necessary
Incorrect temperature display	– The temperature sensor is not immersed deep enough in the measuring solution	– Observe the minimum immersion depth
	– Temperature sensor defective	– Return the D.O. sensor

5 Technical data

General features

Measuring principle	Optical measurement based on photoluminescence
Temperature sensor	Integrated

Dimensions (in mm)

Weight

285 g (with 3 m cable)

Connection cable

Length	3 m
Diameter	5 mm
Smallest allowed bend radius	Fixed installation: 20 mm Flexible use: 60 mm
Plug type	Socket, 4 pins

Pressure resistance

Sensor with connection cable	IP 68 (2.5 x 10 ⁵ Pa or 2.5 bar)
Cable plug	IP 67 (when plugged in)

Measurement conditions

Measuring ranges at 20 °C (68 °F)	0 ... 50 mg/l D.O. 0 ... 500 % D.O. saturation 0 ... 1000 mbar D.O. partial pressure
Temperature range	0 ... 50 °C (32 ... 122 °F)
Immersion depth	min. 65 mm max. 3 m
Operating position	Any
Approach flow	not required

Measuring ranges, accuracy	Variable	Measuring range (max. resolution)	Accuracy (± 1 Digit)
Measuring ranges at 20 °C (68 °F)		0.00 ... 20.00 mg/l O ₂ (± 0.01 mg/l)	± 0.10 mg/l or $\pm 1\%$ of the measured value)
		20.00 ... 50.00 mg/l O ₂ (± 0.01 mg/l)	$\pm 10\%$ of measured value
		0.0 ... 200.0 % O ₂ saturation ($\pm 0.1\%$)	$\pm 1.0\%$ O ₂ saturation or $\pm 1\%$ of the measured value (depending on which value is greater)
		200.0 ... 500.0 % O ₂ saturation ($\pm 1\%$)	$\pm 10\%$ of measured value
		0.0 ... 400.0 mbar O ₂ partial pressure (± 0.1 mbar)	± 2.0 mbar or $\pm 1\%$ of the measured value (depending on which value is greater)
		400.0 ... 1000.0 mbar O ₂ partial pressure (± 0.1 mbar)	$\pm 10\%$ of measured value
Temperature		-5 ... 70 °C (± 0.1 °C)	$\pm 0,2$ °C


Storage conditions

Recommended storing method	In the calibration and storage beaker
Storage temperature	0 ... 50 °C (32 ... 122 °F)
Precision of temperature measurement	± 0.2 K
Working life of the sensor cap	Min. 1 year with authorized use

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VWR International AB
VWR International GmbH
VWR International, LLC