
EN

The replacement Sensor Cap is shipped in a humidified container. The package should not be opened until immediately before Sensor Cap replacement. The Sensor Cap will need to be replaced about once a year, but may last longer. If you are replacing the Sensor Cap, follow the installation and setup instructions below.

- 1 Gasket
- 2 O-Ring
- 3 Sensor Cap

## 1. Sensor Cap replacement

**Caution: Avoid touching the sensing end of the Sensor Cap during the following maintenance procedures. Do not use any tools for these procedures!**

1. Remove the old Sensor Cap assembly by unscrewing it counterclockwise.
2. Carefully remove the o-ring and replace it with the new o-ring included with the replacement Sensor Cap.
3. Apply a thin coat of o-ring lubricant to the clean installed o-ring. Remove any excess o-ring lubricant from the o-ring and/or probe with a lens cleaning tissue.
4. Remove the new Sensor Cap from its humidified container and completely dry the inside cavity of the sensor with a clean lens cleaning tissue. If necessary, clean the clear surface of the sensor on the end of the probe with lens cleaning tissue.
5. Using clockwise motion, thread the Sensor Cap onto the probe assembly until it is finger-tight.
6. After installing the new sensor cap, store the sensor in either water or the grey calibration/storage sleeve with the sponge moistened. If the Sensor Cap dries out, refer to the probe manual for instructions on how to rehydrate the Sensor Cap.

## 2. Enter Sensor Cap coefficients for new Sensor cap

The six coefficients of Calibration Code Label (K1 through K5 and KC) contain the calibration code for this particular Sensor Cap. After installing a new Sensor Cap, connect the probe/cable assembly to the instrument and turn it on. Enter the Sensor Cap coefficients of the Calibration Code Label, included with replacement Sensor Cap, in the instrument menu for *D. O. measurement settings* according to the instrument manual.

If errors are made in entering the Sensor Cap Coefficients, the instrument will block the update. If you continue to get an error message after several entry attempts, contact VWR Technical Support for assistance.

**After entering the Sensor Cap coefficients, perform a 1-point calibration according to instrument manual.**

Important: Be sure to save the Calibration Code label in case you need to reload the calibration coefficients.

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