



Material	C2304-1GAL
Material description	Formic acid
Grade	ACS
CAS Number	64-18-6
Molecular formula	HCO ₂ H
Molecular mass	46.03

Characteristics

Specifications

Assay	Min. 99.0 %
Dilution test	Passes test
SO ₃ (Sulphite)	Passes test
Colour	Max. 15 APHA
Acetic acid	Max. 0.4 %
Heavy metals (as Pb)	Max. 0.001 %
Residue after evaporation	Max. 0.003 %
Cl (Chloride)	Max. 0.001 %
NH ₄ (Ammonium)	Max. 0.005 %
SO ₄ (Sulphate)	Max. 0.003 %
Fe (Iron)	Max. 0.001 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Alicia Bourdeau, Head of the Laboratory - Rouses Point
VWR Chemicals, LLC, 3 Lincoln Blvd, Rouses Point, NY
12979-1045, USA

GF Spec No. MS-001842
GF Part No. CH-PC4103-B1

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.