



Aqua Solutions, Inc.

P.O. Box 70 Deer Park, TX 77536 Ph (281-479-2569)



SPECIFICATION SHEET

SODIUM PHOSPHATE TRIBASIC CRYST / LABORATORY GRADE

Product Code	S3535
Chemical Name / Grade	Sodium Phosphate Tribasic Cryst / Lab
Chemical Formula / CAS#	Na ₃ PO ₄ •12H ₂ O / CAS# 10101-89-0
Appearance	White Crystalline
Harmonization Code	2835.29.3000

CHEMICAL ANALYSIS / SPECIFICATIONS

Assay	90.0% Minimum
Specifications	Minimum / Maximum
P ₂ O ₅ , %	18.5 – 21.0 %
Alkalinity weight, %	18.0 – 20.5 %
Loss On Ignition	45.00 – 57.00 %

This product is suitable for use in applications where the grade indicated above is required.

PREPARATION DATE (MM/DD/YY)

Original 0.0	04/13/2018: STN
Revision	

PLEASE REFER TO THE MATERIAL SAFETY DATA SHEET FOR HEALTH, SAFETY, AND REGULATORY INFORMATION ON THIS PRODUCT.