



Aqua Solutions, Inc.

P.O. Box 70 Deer Park, TX 77536 Ph (281-479-2569)



SPECIFICATION SHEET

SODIUM ACETATE / REAGENT ACS

Product Code	S2420
Chemical Name / Grade	Sodium Acetate / Reagent ACS
Chemical Formula / CAS#	CH ₂ COONa / CAS# 127-09-3
Appearance	White deliquescent powder
Country of Origin	China
Harmonization Code	2915.29.2000

CHEMICAL ANALYSIS / SPECIFICATIONS

Assay	99.0% min
Specifications	Maximum
Loss on Drying	1.0%
Heavy Metals (as Pb)	0.001%
Chloride (Cl)	0.002%
Sulfate (SO ₄)	0.003%
Calcium (Ca)	0.005%
Insoluble matter	0.01%
Phosphate (PO ₄)	0.001%
Magnesium (Mg)	0.002%
Iron (Fe)	0.001%

This product is suitable for use in applications where the grade indicated above is required.

PREPARATION DATE (MM/DD/YY)

Original 0.0	06/16/2014: STN
Revision	

PLEASE REFER TO THE MATERIAL SAFETY DATA SHEET FOR HEALTH, SAFETY, AND REGULATORY INFORMATION ON THIS PRODUCT.