



Aqua Solutions, Inc.

P.O. Box 70 Deer Park, TX 77536 Ph (281-479-2569)



SPECIFICATION SHEET

BARIUM CHLORIDE 20% W/V SOLUTION

Product Code	0920
Product Name	Barium Chloride 20% w/v Solution
Actual Assay	20.0 +/- 0.1% Barium Chloride
Prepared to Formulation on File	Yes
Appearance and Color	Clear-Slightly Hazy Liquid
Country of Origin	USA
Harmonization Code	2827.39.4500

The weights and /or volumes used to prepare this product are N.I.S.T. Traceable. All balances used in the preparation of product are calibrated daily against N.I.S.T. Traceable weights. The balances are maintained and serviced on a regular basis by an outside certified company. All Volumetric glassware used is N.I.S.T. Traceable and certified as meeting Class A specifications.

PREPARATION DATE (MM/DD/YY)

Original 0.0	04/12/2018 STN
Revision	

PLEASE REFER TO THE MATERIAL SAFETY DATA SHEET FOR HEALTH, SAFETY, AND REGULATORY INFORMATION ON THIS PRODUCT.